




Application of EUROIMMUN tests for COVID-19 diagnostics



Correct application of our tests

- Identification of an acute infection by direct pathogen detection using PCR and antigen ELISA: **EURORealTime SARS-CoV-2, EURORealTime SARS-CoV-2/Influenza A/B** as well as **SARS-CoV-2 Antigen ELISA**
- Further monitoring of patients with acute infection (confirmed by positive PCR): Antibody detection using the **Anti-SARS-CoV-2 ELISA (IgA)** and **Anti-SARS-CoV-2 NCP ELISA (IgM)**
- Detection of a past contact with the pathogen (from 2 to 4 weeks after infection): Antibody detection using the **Anti-SARS-CoV-2 ELISA (IgG)** or **Anti-SARS-CoV-2 QuantiVac ELISA (IgG)** and **Anti-SARS-CoV-2 NCP ELISA (IgG)**



Areas of application of EUROIMMUN SARS-CoV-2 tests

	PCR	Antigen ELISA	IgA ELISA	IgM ELISA	IgG ELISA
Detection of acute infections	✓	✓	✗	✗	✗
Screening for detection of acute infections	✓	✓	✗	✗	✗
Detection of a pathogen contact (up to day 10)	✓	✓	✗/✓	✗/✓	✗
Detection of a pathogen contact (from week 2 to 4)	✗	✗	✗	✗	✓
Further monitoring after acute infection (confirmed by direct detection)	✓	✓	✓	✓	✓
Detection of past infections	✗	✗	✗	✗	✓

Interpretation of SARS-CoV-2 test results

	PCR	Antigen ELISA	IgA ELISA	IgM ELISA	IgG ELISA
Early phase of acute infection (up to day 10)	+	+	+/-	+/-	- (+)
Late phase of acute infection (from day 10)	- (+)	-	+/-	+/-	+ (-)
Past infection	-	-	+/-	-	+
No statement possible (follow-up sample recommended)	-	-	+	+	-

PCR: EURORealTime SARS-CoV-2, EURORealTime SARS-CoV-2/Influenza A/B; antigen ELISA: SARS-CoV-2 Antigen ELISA; IgA ELISA: Anti-SARS-CoV-2 ELISA (IgA); IgM ELISA: Anti-SARS-CoV-2 NCP ELISA (IgM); IgG ELISA: Anti-SARS-CoV-2 ELISA (IgG), Anti-SARS-CoV-2 QuantiVac ELISA (IgG), Anti-SARS-CoV-2 NCP ELISA (IgG)

Characteristics of EUROIMMUN SARS-CoV-2 tests

- **EURORealTime SARS-CoV-2:** High sensitivity and reliability due to simultaneous detection of two target sequences of SARS-CoV-2; only one reaction per sample
- **EURORealTime SARS-CoV-2/Influenza A/B:** For direct detection of SARS-CoV-2 and influenza viruses (types A and B); for differential diagnostic clarification of symptoms that can be associated with COVID-19 as well as influenza
- **Anti-SARS-CoV-2 ELISA (IgG):** High specificity of 99.6% due to the use of the spike protein domain S1 including the immunologically relevant receptor binding domain (RBD) – the main target antigen for virus-neutralising antibodies
- **Anti-SARS-CoV-2 QuantiVac ELISA (IgG):** Quantitative determination of Anti-SARS-CoV-2 S1/RBD antibodies in standardised units based on the “First WHO International Standard for anti-SARS-CoV-2 immunoglobulin” (NIBSC code: 20/136); excellent correlation with neutralisation tests; supports the detection and assessment of the immune response following infection or vaccination
- **Anti-SARS-CoV-2 NCP ELISA (IgM or IgG):** High specificity of the ELISAs of 98.6% (IgM) and 99.8% (IgG), due to the use of a modified nucleocapsid protein (NCP) containing exclusively diagnostically relevant epitopes
- **Anti-SARS-CoV-2 ELISA (IgA):** Use of the S1 domain of the spike protein as the antigen; for supportive analysis of follow-up samples after infection confirmed by PCR; indicates a beginning immune reaction