




# Antigen assay for acute diagnostics of COVID-19 in the laboratory



## SARS-CoV-2 Antigen ELISA

- Laboratory diagnostic assay for direct pathogen detection by semiquantitative determination of the SARS-CoV-2-specific nucleocapsid protein in swab samples of the upper respiratory tract
- Supporting acute diagnostics of SARS-CoV-2 infections – mainly during COVID-19 outbreaks
- High correlation with PCR test results: 100% negative agreement (specificity) and 93.6% positive agreement (sensitivity)
- Established ELISA method – suitable for each diagnostic laboratory and completely automatable
- Extension of the test application for saliva samples currently under validation

**Order number: EQ 2606-9601**