



A powerful trio from EUROIMMUN

Complete analysis of the immune response after SARS-CoV-2 infection or vaccination



Anti-SARS-CoV-2 QuantiVac (IgG)

ELISA for quantification of the anti-S1-/RBD IgG antibody concentration in standardised units (BAU/ml)



SARS-CoV-2 NeutraLISA

Neutralisation test in ELISA format for detection of an inhibitory effect of neutralising antibodies directed against S1/RBD



SARS-CoV-2 IGRA*

ELISA-based interferon gamma release assay (IGRA) for determination of the activity of SARS-CoV-2-reactive T cells by stimulation through S1/RBD

S1: S1 domain of the spike protein of SARS-CoV-2; **RBD:** receptor binding domain

* Already available for research use (not for in vitro diagnostics in the sense of EU Directive 98/79/EC)
soon also available with CE mark